

320ml / 10oz (80mm) White Single Wall BioCup

Product: 320ml / 10oz (80mm) White Single Wall BioCup

SKU: BC-10W

URL: <https://www.sipandpak.com.au/product/320ml-10oz-80mm-white-single-wall-biocup/>

BioPak specification

Unit weight (g)	9.2
Packaging type	Cup
Raw material	Paper - FSC(TM) Mix certified paper - Bioplastic - Ingeo®, PLA lining/coating
Lining	Bioplastic - Ingeo®, PLA
Home compostable	No
Industrial compostable	Certified to AS4736 Australian Standards, Certified to EN13432 European Standard
Recyclable	Not Recyclable
Certifications	Industrially compostable: Certified to AS4736 Australian Standards, Certified to EN13432 European Standard · Environmental production certified: ISO 14001 · Factory food safety certified: BRC · FSC (material)
PFAS	Not applicable

Impact estimate (500 units/week x 52 weeks)

Annual units	26,000
Format benchmark	12oz single-wall cup (conventional PE-lined benchmark)
Industry benchmark	12oz single-wall cup (conventional PE-lined benchmark)
Sustainable unit mass (g)	9.2
Conventional unit mass (g)	13.0
Fossil polymer saved (kg/year)	26.00
Mass saved (kg/year)	98.80
Polymer reduction (%)	100.0
Mass reduction (%)	29.2
Conventional packaging total (kg/year)	338.00
Sustainable packaging total (kg/year)	239.20
Confidence	High
Data source	sku
Sustainability highlights	Paper - FSC(TM) Mix certified paper - Bioplastic - Ingeo®, PLA lining/coating; No; Industrially compostable: Certified to AS4736 Australian Standards, Certified to EN13432 European Standard · Environmental production certified: ISO 14001 · Factory food safety certified: BRC · FSC (material); Recyclable: Not Recyclable

Methodology

Equivalent Format Comparison (EFC): BioPak specification unit mass vs industry conventional packaging mass for the matched format, plus fossil polymer estimate using typical PE/PS/PP fractions. SKU BC-10W weight from BioPak specifications. Indicative estimate only - not a certified life-cycle assessment (ISO 14040/14044). Compares BioPak published unit masses to conventional equivalent-format benchmarks used in Australian foodservice (PE-lined cups, PP clamshells, PS lids). For procurement reporting, request substantiated SKU documentation from Sip and Pak. Typically within $\pm 15\%$ for fossil polymer estimate when a BioPak SKU and format benchmark both match.