

390ml / 12oz (90mm) Christmas Single Wall BioCup

Product: 390ml / 12oz (90mm) Christmas Single Wall BioCup

SKU: BC-12-CHRISTMAS

URL: <https://www.sipandpak.com.au/product/390ml-12oz-90mm-christmas-single-wall-biocup/>

BioPak specification

Unit weight (g)	9.0
Packaging type	Cup
Raw material	Non-FSC Paper - Bioplastic - PLA lining/coating
Lining	Bioplastic - PLA
Home compostable	No
Industrial compostable	Certified to AS4736 Australian Standards, Certified to EN13432 European Standard
Recyclable	Not Recyclable
Certifications	Industrially compostable: Certified to AS4736 Australian Standards, Certified to EN13432 European Standard · Environmental production certified: ISO 14001 · Factory food safety certified: FSSC 22000, ISO 22000
PFAS	Not applicable

Impact estimate (500 units/week x 52 weeks)

Annual units	26,000
Format benchmark	12oz single-wall cup (conventional PE-lined benchmark)
Industry benchmark	12oz single-wall cup (conventional PE-lined benchmark)
Sustainable unit mass (g)	9.0
Conventional unit mass (g)	13.0
Fossil polymer saved (kg/year)	26.00
Mass saved (kg/year)	104.00
Polymer reduction (%)	100.0
Mass reduction (%)	30.8
Conventional packaging total (kg/year)	338.00
Sustainable packaging total (kg/year)	234.00
Confidence	High
Data source	sku
Sustainability highlights	Non-FSC Paper - Bioplastic - PLA lining/coating; No; Industrially compostable: Certified to AS4736 Australian Standards, Certified to EN13432 European Standard · Environmental production certified: ISO 14001 · Factory food safety certified: FSSC 22000, ISO 22000; Recyclable: Not Recyclable

Methodology

Equivalent Format Comparison (EFC): BioPak specification unit mass vs industry conventional packaging mass for the matched format, plus fossil polymer estimate using typical PE/PS/PP fractions. SKU BC-12-CHRISTMAS weight from BioPak specifications. Indicative estimate only - not a certified life-cycle assessment (ISO 14040/14044). Compares BioPak published unit masses to conventional equivalent-format benchmarks used in Australian foodservice (PE-lined cups, PP clamshells, PS lids). For procurement reporting, request substantiated SKU documentation from Sip and Pak. Typically within $\pm 15\%$ for fossil polymer estimate when a BioPak SKU and format benchmark both match.