

740ml / 24oz Kraft BioBowl

Product: 740ml / 24oz Kraft BioBowl

SKU: BSCK-24-GS-V2

URL: <https://www.sipandpak.com.au/product/740ml-24oz-kraft-biobowl/>

BioPak specification

Unit weight (g)	16.7
Packaging type	Container
Raw material	Paper - FSC(TM) Mix certified paper - Bioplastic - PLA lining/coating
Lining	Bioplastic - PLA
Home compostable	No
Industrial compostable	Certified to AS4736 Australian Standards, Certified to EN13432 European Standards
Recyclable	Not Recyclable
Certifications	Industrially compostable: Certified to AS4736 Australian Standards, Certified to EN13432 European Standards · Environmental production certified: ISO 14001 · Factory food safety certified: BRC · FSC (material)
PFAS	Not applicable

Impact estimate (500 units/week x 52 weeks)

Annual units	26,000
Format benchmark	Food bowl (conventional PP/paper benchmark)
Industry benchmark	Food bowl (conventional PP/paper benchmark)
Sustainable unit mass (g)	16.7
Conventional unit mass (g)	18.0
Fossil polymer saved (kg/year)	416.00
Mass saved (kg/year)	33.54
Polymer reduction (%)	100.0
Mass reduction (%)	7.2
Conventional packaging total (kg/year)	468.00
Sustainable packaging total (kg/year)	434.46
Confidence	High
Data source	sku
Sustainability highlights	Paper - FSC(TM) Mix certified paper - Bioplastic - PLA lining/coating; No; Industrially compostable: Certified to AS4736 Australian Standards, Certified to EN13432 European Standards · Environmental production certified: ISO 14001 · Factory food safety certified: BRC · FSC (material); Recyclable: Not Recyclable

Methodology

Equivalent Format Comparison (EFC): BioPak specification unit mass vs industry conventional packaging mass for the matched format, plus fossil polymer estimate using typical PE/PS/PP fractions. SKU BSCK-24-GS-V2 weight from BioPak specifications. Indicative estimate only - not a certified life-cycle assessment (ISO 14040/14044). Compares BioPak published unit masses to conventional equivalent-format benchmarks used in Australian foodservice (PE-lined cups, PP clamshells, PS lids). For procurement reporting, request substantiated SKU documentation from Sip and Pak. Typically within $\pm 15\%$ for fossil polymer estimate when a BioPak SKU and format benchmark both match.